



WHILL Model F

The perfect travel chair has arrived.



The Perfect Travel Power Chair Has Arrived.

Lightweight, foldable, maneuverable. The Model F allows you to explore your world with greater mobility and independence.

Experience the same performance and smart technology that WHILL is known for with a new foldable design.



Foldable Design

Easily transport the Model F on buses, trains, and planes – or in the trunk of a car or taxi.*



When not in use, the Model F can be conveniently stored in a closet or out-of-the-way corner of the room.



*Model F is not to be used as a seat while it is being transported.

The information contained herein is correct at the time of printing, and we reserve the right to alter specifications without prior notice.

Lightweight

Weighing only 52.9 lbs., the Model F is one of the lightest power chairs on the market today.**



Top Speed



Range



Obstacle Clearance



30.7" turning radius for tight turns

Connect Intelligently

The WHILL app allows you to check key device information and drive the Model F remotely.



Drive remotely, lock/unlock, and manage speed settings via the WHILL smartphone app.

Drive with Peace of Mind

Stay engaged in activities and embrace your surroundings knowing the Model F will keep you safe while you ride.



Safely drives up or down 10° slopes.



**Weight is 58.9 lbs. if battery is not removed. Battery weighs 6 lbs.

Note: Driving range of 12.4 miles will vary due to user weight, drive surface, terrain, ambient temperature, speed setting, and battery charge and condition.

Customization

Choose from five colors and a great selection of accessories to match your personal style and enhance your driving experience.

Color Options



Controller can be positioned on either the left or right side.

Accessories

◦ Non-medical options



Smart Key



Joystick



Cane Holder



Waterproof Cover◦



Side Bag◦



Controller Cover◦



Clothes Guard◦



Heel Stopper◦

Specifications

Driving Range	Up to 12.4 Miles	Maximum Speed***	3.7 MPH
Foldable	Yes	Max Incline	10°
Turning Radius	30.7"	Obstacle Clearance	1.4"
Seat Sizes Available	18"	Width Between Armrests	Adjustable: 2 positions (16", 18")
Device Width	21.9" / 23.8"	Device Length	36.8" (18.3" when folded)
Device Weight	58.9 lbs	Charging Time	~ 5 hours
App Enabled	Yes, with remote control	Air Travel Approved	Yes
Max Weight Capacity	Up to 250 lbs	Swappable Arm Covers	5 Colors
Battery	Lithium-Ion (25.3V 10Ah)	Water Resistance	IPX4

***Acceleration and/or deceleration may change when driving on steep slopes.

WHILL Inc.

303 Twin Dolphin Dr., 6th Floor, Redwood City, CA 94065
 (844) 699-4455 us.info@whill.inc www.whill.inc



Indications for Use: The intended use of the WHILL Model F powered wheelchair is to provide outdoor and indoor mobility to persons limited to a seated position that are capable of operating a powered wheelchair.

Wheelchair Class: A FDA Compliance: The WHILL Model F has received 510(k) clearance for marketing as a Class II medical device.

740-00236