

WHILL Model Fi



WHILL

The Perfect Travel Power Chair Has Arrived.

Lightweight, foldable, maneuverable. The Model Fi allows you to explore your world with greater mobility and independence.

Experience the same performance and smart technology that WHILL is known for with a new foldable design.



Foldable Design

Easily transport the Model Fi on buses, trains, and planes – or in the trunk of a car or taxi.*



When not in use, the Model Fi can be conveniently stored in a closet or out-of-the-way corner of the room.



*Model Fi is not to be used as a seat while it is being transported.

The information contained herein is correct at the time of printing, and we reserve the right to alter specifications without prior notice.

Lightweight

Weighing only 52.9 lbs., the Model Fi is one of the lightest power chairs on the market today.**



Top Speed



Range



Obstacle Clearance



30.7" turning radius for tight turns

Connect Intelligently

The WHILL app allows you to check key device information and drive the Model Fi remotely.



Drive remotely, lock/unlock, and manage speed settings via the WHILL smartphone app.

Drive with Peace of Mind

Stay engaged in activities and embrace your surroundings knowing the Model Fi will keep you safe while you ride.



Safely drives up or down 10° slopes.



**Weight is 58.9 lbs. if battery is not removed. Battery weighs 6 lbs.

Note: Driving range of 12.4 miles will vary due to user weight, drive surface, terrain, ambient temperature, speed setting, and battery charge and condition.

Customization

Choose from five colors and a great selection of accessories to match your personal style and enhance your driving experience.

Color Options



- White ● Black ● Light Green
- Red ● Light Blue



Controller can be positioned on either the left or right side.

Accessories



Battery



Battery
Charger



Lap Belt



Smart Key



Joystick



Basket



Cane Holder



Dust Cover



Side Bag

Specifications

Driving Range	Up to 12.4 Miles
Foldable	Yes
Turning Radius	30.7"
Seat Sizes Available	18"
Device Width	21.9" / 23.8"
Device Weight	58.9 lbs
App Enabled	Yes, with remote control
Max Weight Capacity	Up to 250 lbs
Battery	Lithium-Ion [25.3V 10Ah]

Maximum Speed***	3.7 MPH
Max Incline	10°
Obstacle Clearance	1.4"
Width Between Armrests	Adjustable: 2 positions (16", 18")
Device Length	36.8" (18.3" when folded)
Charging Time	~ 5 hours
Air Travel Approved	Yes
Swappable Arm Covers	5 Colors
Water Resistance	IPX4

***Acceleration and/or deceleration may change when driving on steep slopes.

WHILL Inc.

951 Mariners Island Blvd., Suite 300, San Mateo, CA 94404
(844) 699-4455 us.info@whill.inc www.whill.inc

