



Product Comparison



Model Ci



Model A

MSRP/MAP	\$4,499 / \$3,999	\$9,999 / \$7,999
Product Versatility		
Drive System	2WD	4WD
Number of Motors	2	2
Incline/Decline Capability	10 degrees/17.6%	10 degrees/17.6%
Obstacle Clearance (e.g. curbs, rocks)	2 inches	3 inches
Front Wheels	9.8" diameter omni-wheels	9.8" diameter omni-wheels
Rear Wheels	10.5" cushioned solid tires	12.5" cushioned solid tires
Maximum Speed	5.0 mph	5.5 mph
Range on a Full Charge	Up to 10 miles	Up to 12 miles
Battery Type and Quantity	One, LI Ion	Two, sealed glass mat
Swappable Batteries	✓	✗
Turning Radius	29.9"	28"
Driving, Steering, Maneuverability		
Controller Types	Mouse and Joystick	Mouse and Easy Grip options
Number of Speed Settings	4	Fast, medium, slow
Controller Location	RT or LT, user changeable	RT or LT, performed by reseller
iPhone App for Remote Driving	✓	✓
iPhone and Android App to Modify Speed Settings	✓	✓
Range Indicator on Control Pad	Displayed 1-100%	1-5 bars
Braking System	Electromagnetic	Electromagnetic
Integrated Horn	✓	✗
User Specifications		
18" Seat Width + Cushion Depth Options	✓ 18" cushion depth	✗
16" Seat Width + Cushion Depth Options	✓ with 16", 18" depth options	✓ 16", 17" 18", 19", 20" seat depth
Total weight limit, including any carried items	250LB	220LB
Lap Restraint	Basic	Basic
Storage Under Seat	✓ 5 gallons	✗
Storage Under Armrests	Custom sidebag available	✗
Attach Accessories (cupholder, iphone, etc.)	✓	✓
Powered Seat Slide Up to 5.9" for Extended Reach at Tables and Desks	✗	✓
User Height Specifications	From 4'8"-6'2"	From 4'-6'2"
Other Product Specifications		
USB Port	1-USB port	✗
Color Options	8 color options	White or black
Armrests	✓	✗
Product Width	16" Seat = 21.5" / 18" Seat = 23.5"	16" seat/23.6", 18" seat/25.7"
Product Length	38.8"	35" - 39"
Product Height (without backrest cushion on)	29"	31" - 33"
Product Weight	16" Seat = 112 LB / 18" Seat = 114 LB	250lbs
Wheelbase	23.5"	20.6"
FDA-Cleared	✗	✗
Rehab Benefits		
Seat Cushion	General use - 16", 18" [deep foam]	General use - 16", 18", 20" [deep foam]
Back Support Options	1-size, foam filled	Low, mid foam filled
Pressure Relief	✗	✗
Back Post Expansions	✗	✗
4-Point Tie-Down Option	✗	✗
Q'Straint Compatible	✗	✗
Partially Reimbursable by Insurance	✗	✗